

Anti-Calcyclin CACYBP Monoclonal Antibody

Catalog Number: M04327


Overview

Product Name	Anti-Calcyclin CACYBP Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Calcyclin CACYBP Monoclonal Antibody catalog # M04327. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 7D11
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9HB71


Technical Details

Immunogen	Purified recombinant fragment of Calcyclin expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-Calcyclin CACYBP Monoclonal Antibody (M04327) Images



Western blotting validation for Anti-Calcyclin CACYBP Monoclonal Antibody M04327 Western Blot (WB) analysis using Calcyclin Monoclonal Antibody against truncated calcyclin recombinant protein. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Immunohistochemistry validation of CACYBP using Anti-Calcyclin CACYBP Monoclonal Antibody (M04327). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Brain glioma tissue showing nuclear localization with DAB staining using Calcyclin Monoclonal Antibody. For more protocol information of IHC

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calcyclin CACYBP Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.