

Anti-TNFRSF21 / DR6 / CD358 Reference Antibody (Abbvie patent anti-TNFRSF21)

Catalog Number: M04348-2


Overview

Product Name	Anti-TNFRSF21 / DR6 / CD358 Reference Antibody (Abbvie patent anti-TNFRSF21)
Reactive Species	Human
Description	Boster Bio Anti-TNFRSF21 / DR6 / CD358 Reference Antibody (Abbvie patent anti-TNFRSF21) (Catalog # M04348-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.857EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	O75509


Technical Details

Isotype	IgG1
Form	Liquid

Anti-TNFRSF21 / DR6 / CD358 Reference Antibody (Abbvie patent anti-TNFRSF21) (M04348-2) Images



Anti-TNFRSF21/DR6/CD358 Reference Antibody (Abbvie patent anti-TNFRSF21) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-TNFRSF21/DR6/CD358 Reference Antibody (Abbvie patent anti-TNFRSF21) is more than 95%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFRSF21 / DR6 / CD358 Reference Antibody (Abbvie patent anti-TNFRSF21)

For Research Use Only. Not for use in diagnostic procedures.