

Anti-NCOA4 Antibody (Monoclonal, 35N05)

Catalog Number: M04368-3

About NCOA4

Nuclear receptor coactivator 4 is a protein that in humans is encoded by the NCOA4 gene. This gene encodes an androgen receptor coactivator. The encoded protein interacts with the androgen receptor in a ligand-dependent manner to enhance its transcriptional activity. Chromosomal translocations between this gene and the ret tyrosine kinase gene, also located on chromosome 10, have been associated with papillary thyroid carcinoma. Alternatively spliced transcript variants have been described. Pseudogenes are present on chromosomes 4, 5, 10, and 14.

Overview

Product Name	Anti-NCOA4 Antibody (Monoclonal, 35N05)
Reactive Species	Human
Description	Boster Bio Anti-NCOA4 Antibody (Monoclonal, 35N05) catalog # M04368-3. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 35N05
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13772

Technical Details

Immunogen	Recombinant protein within human NCOA4 aa 3-398.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NCOA4 Antibody (Monoclonal, 35N05)

For Research Use Only. Not for use in diagnostic procedures.