

Anti-Oct-2 POU2F2 Monoclonal Antibody

Catalog Number: M04407-1

Overview

Product Name	Anti-Oct-2 POU2F2 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Oct-2 POU2F2 Monoclonal Antibody catalog # M04407-1. Tested in IHC applications. This antibody reacts with Human.
Application	IHC
Clonality	Monoclonal AS1A36
Formulation	Liquid form in PBS containing 50% glycerol, 0.5% BSA & 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P09086

Technical Details

Immunogen	Synthesized peptide derived from human Oct-2 (Oct2)
Form	Liquid
Concentration	1 mg/ml
Purification	Oct-2 (Oct2) (AS1A36) antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used IHC 1:500</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Oct-2 POU2F2 Monoclonal Antibody