

Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378

Catalog Number: M04431-1


Overview

Product Name	Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378
Reactive Species	Human
Description	Boster Bio Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378 (Catalog # M04431-1). Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal RM378
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16633


Technical Details

Immunogen	A peptide corresponding to the C-terminus of human BOB-1/OBF-1
Cross Reactivity	This antibody reacts to human BOB-1/OBF-1 (B-cell-specific coactivator) (POU domain class 2-associating factor 1)
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	Immunohistochemistry (IHC): 1:1000-1:2000 dilution WB: 1:5000-1:10000 dilution.

Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378 (M04431-1) Images



Western Blotting result Western Blot of Daudi cells lysate using anti-BOB-1/OBF-1 rabbit monoclonal antibody (Clone RM378) at a 1:10000 dilution.



IHC result Immunohistochemical staining of formalin fixed and paraffin embedded human tonsil tissue section using anti-BOB-1/OBF-1 rabbit monoclonal antibody (Clone RM378) at a 1:1250 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378

For Research Use Only. Not for use in diagnostic procedures.