

Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378

Catalog Number: M04431-1

Overview

Product Name	Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378
Reactive Species	Human
Description	Boster Bio Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378 (Catalog # M04431-1). Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal RM378
Formulation	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16633

Technical Details

Immunogen	A peptide corresponding to the C-terminus of human BOB-1/OBF-1
Cross Reactivity	This antibody reacts to human BOB-1/OBF-1 (B-cell-specific coactivator) (POU domain class 2-associating factor 1)
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>Immunohistochemistry (IHC): 1:1000-1:2000 dilution</p> <p>WB: 1:5000-1:10000 dilution.</p>

Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378 (M04431-1) Images




Figure 1. Western Blotting result
Western Blot of Daudi cells lysate using anti-BOB-1/OBF-1 rabbit monoclonal antibody (Clone RM378) at a 1:10000 dilution.




Figure 2. IHC result
Immunohistochemical staining of formalin fixed and paraffin embedded human tonsil tissue section using anti-BOB-1/OBF-1 rabbit monoclonal antibody (Clone RM378) at a 1:1250 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BOB-1/OBF-1 Rabbit Monoclonal Antibody, Clone#RM378