

Anti-SERPINB1 Mouse Monoclonal Antibody [Clone ID: OTI1D1]

Catalog Number: M04491-1


Overview

Product Name	Anti-SERPINB1 Mouse Monoclonal Antibody [Clone ID: OTI1D1]
Reactive Species	Dog, Human, Monkey
Description	Boster Bio SERPINB1 mouse monoclonal antibody, clone OTI1D1 (formerly 1D1). Catalog# M04491-1. Tested in FC, IF, WB. This antibody reacts with Human, Monkey, Dog.
Application	Flow Cytometry, IF, WB
Clonality	Monoclonal OTI1D1
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P30740


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 77-306 of human SERPINB1 (NP_109591) produced in E.coli.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500~2000 IF 1:100 Flow Cytometry 1:100


Anti-SERPINB1 Mouse Monoclonal Antibody [Clone ID: OTI1D1] (M04491-1) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SERPINB1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SERPINB1.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-SERPINB1 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



HEK293T cells transfected with either SERPINB1 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-SERPINB1 antibody (M04491-1)



Flow cytometric Analysis of Jurkat cells



Anti-SERPINB1 mouse monoclonal antibody (M04491-1) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY SERPINB1.

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SERPINB1 Mouse Monoclonal Antibody [Clone ID: OT11D1]

For Research Use Only. Not for use in diagnostic procedures.