

Anti-Forkhead box protein D3 FoxD3 Monoclonal Antibody

Catalog Number: M04495-1


Overview

Product Name	Anti-Forkhead box protein D3 FoxD3 Monoclonal Antibody
Reactive Species	Human, Monkey, Mouse
Description	Boster Bio Anti-Forkhead box protein D3 FoxD3 Monoclonal Antibody catalog # M04495-1. Tested in ELISA, IHC, WB applications. This antibody reacts with Human, Monkey, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 5G9
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9UJU5


Technical Details

Immunogen	Purified recombinant fragment of human FoxD3 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-Forkhead box protein D3 FoxD3 Monoclonal Antibody (M04495-1) Images



Western Blot analysis using FoxD3 Monoclonal Antibody against NTERA-2 (1), HUVE-12 (2), HEK293 (3), HeLa (4), Jurkat (5), K562 (6), RAW264.7 (7), NIH/3T3 (8), and COS7 (9) cell lysate.



Immunohistochemistry analysis of paraffin-embedded lung cancer tissues (left) and ovarian cancer tissues (right) with DAB staining using FoxD3 Monoclonal Antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Forkhead box protein D3 FoxD3 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.