

Anti-Glycine decarboxylase/GLDC Antibody Picoband® (monoclonal, 3D3D3) Biotin Conjugated

Catalog Number: M04777-Biotin

About GLDC

Glycine decarboxylase also known as glycine cleavage system P protein or glycine dehydrogenase is an enzyme that in humans is encoded by the GLDC gene. Degradation of glycine is brought about by the glycine cleavage system, which is composed of four mitochondrial protein components: P protein (a pyridoxal phosphate-dependent glycine decarboxylase), H protein (a lipoic acid-containing protein), T protein (a tetrahydrofolate-requiring enzyme), and L protein (a lipoamide dehydrogenase). The protein encoded by this gene is the P protein, which binds to glycine and enables the methylamine group from glycine to be transferred to the T protein. Defects in this gene are a cause of nonketotic hyperglycinemia (NKH).

Overview

| | |
|----------------------|--|
| Product Name | Anti-Glycine decarboxylase/GLDC Antibody Picoband® (monoclonal, 3D3D3) Biotin Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | WB, IHC, ELISA |
| Clonality | Monoclonal 3D3D3 |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Mouse |
| Uniprot ID | P23378 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human Glycine decarboxylase/GLDC recombinant protein (Position: K574-S1020). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Mouse IgG1 |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Biotin |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by |

end users.
ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Glycine decarboxylase/GLDC Antibody (monoclonal, 3D3D3) - Biotin

For Research Use Only. Not for use in diagnostic procedures.