

Anti-Hyaluronan synthase 1 HAS1 Monoclonal Antibody

Catalog Number: M04784


Overview

Product Name	Anti-Hyaluronan synthase 1 HAS1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Hyaluronan synthase 1 HAS1 Monoclonal Antibody catalog # M04784. Tested in ELISA, IF, WB applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Monoclonal 3E10
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q92839


Technical Details

Immunogen	Purified recombinant fragment of human HAS1 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000

Anti-Hyaluronan synthase 1 HAS1 Monoclonal Antibody (M04784) Images



Western Blot analysis using HAS1 Monoclonal Antibody against recombinant protein of human HAS1 (aa70-243).



Immunofluorescence analysis of U251 cells using HAS1 Monoclonal Antibody (green). Red: Actin filaments have been labeled with DY-554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Hyaluronan synthase 1 HAS1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.