

## Anti-Annexin A3 Antibody Picoband® (monoclonal, 2H3H8) Fluoro647 Conjugated

Catalog Number: M04796-2-Fluoro647

### About ANXA3

Annexin A3 is a protein that in humans is encoded by the Annexin A3 gene. The Annexin A3 gene contains 13 exons and spans 58 kb of genomic DNA. The Annexin A3 gene is mapped to 4q21. It is abnormally expressed in fetuses of both IVF and ICSI, which may contribute to the increase risk of birth defects in these ART. This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions in the inhibition of phospholipase A2 and cleavage of inositol 1,2-cyclic phosphate to form inositol 1-phosphate. This protein may also play a role in anti-coagulation.

### Overview

Product Name	Anti-Annexin A3 Antibody Picoband® (monoclonal, 2H3H8) Fluoro647 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Monoclonal 2H3H8
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	P12429

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human Annexin A3, different from the related mouse sequence by one amino acid, and from the related rat sequence by three amino acids.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2b
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Annexin A3 Antibody (monoclonal, 2H3H8) - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.