

## Anti-LRRFIP1 monoclonal Antibody

Catalog Number: M04854


### Overview

Product Name	Anti-LRRFIP1 monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-LRRFIP1 monoclonal Antibody catalog # M04854. Tested in WB applications. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Monoclonal 3E10
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q32MZ4

### Technical Details

Immunogen	Recombinant full length Human LRRFIP1.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000

## Anti-LRRFIP1 monoclonal Antibody (M04854) Images



Western blot detection of LRRFIP1 in Raji, C2C12, MCF7&Jurkat cell lysates using LRRFIP1 antibody (1:1000 diluted) .

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LRRFIP1 monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.