

Anti-PKA 2 beta PRKAR2B Rabbit Monoclonal Antibody

Catalog Number: M04943

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-PKA 2 beta PRKAR2B Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PKA 2 beta PRKAR2B Rabbit Monoclonal Antibody catalog # M04943. Tested in WB, IHC, ICC/IF, IP, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Monoclonal ABBB-16
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P31323


Technical Details

Immunogen	A synthesized peptide derived from human PKA 2 beta
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography

Suggested Dilutions

WB 1:500-2000
IHC 1:50-200
ICC/IF 1:50-200
IP 1:20
FC 1:20

Anti-PKA 2 beta PRKAR2B Rabbit Monoclonal Antibody (M04943) Images



Western blot analysis of PKA 2 beta expression in human fetal brain lysate.

[Anti-PKA 2 beta PRKAR2B Rabbit Monoclonal Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.