

Anti-RGMA Reference Antibody (elezanumab)

Catalog Number: M04984-3


Overview

Product Name	Anti-RGMA Reference Antibody (elezanumab)
Reactive Species	Cynomolgus monkey, Human, Mouse, Rat
Description	Boster Bio Anti-RGMA Reference Antibody (elezanumab) (Catalog # M04984-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Cynomolgus monkey, Human, Mouse. Endotoxin: < 0.872EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q96B86


Technical Details

Isotype	IgG1
Form	Liquid

Anti-RGMA Reference Antibody (elezanumab) (M04984-3) Images



Anti-RGMA Reference Antibody (elezanumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-RGMA Reference Antibody (elezanumab) is more than 100%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RGMA Reference Antibody (elezanumab)

For Research Use Only. Not for use in diagnostic procedures.