

## Anti-Fatty aldehyde dehydrogenase ALDH3A2 Monoclonal Antibody

Catalog Number: M05079-1


### Overview

Product Name	Anti-Fatty aldehyde dehydrogenase ALDH3A2 Monoclonal Antibody
Reactive Species	Human, Mouse, Rabbit, Rat
Description	Boster Bio Anti-Fatty aldehyde dehydrogenase ALDH3A2 Monoclonal Antibody catalog # M05079-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rabbit, Rat.
Application	WB
Clonality	Monoclonal 9E9
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P51648

### Technical Details

Immunogen	Purified recombinant human ALDH3A2 protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG1
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000

## Anti-Fatty aldehyde dehydrogenase ALDH3A2 Monoclonal Antibody (M05079-1) Images



Western Blot analysis using ALDH3A2 Monoclonal Antibody against A549, U251, HepG2 cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fatty aldehyde dehydrogenase ALDH3A2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.