

## Anti-OXSR1 Mouse Monoclonal Antibody [Clone ID: OTI1F3]

Catalog Number: M05141


### Overview

Product Name	Anti-OXSR1 Mouse Monoclonal Antibody [Clone ID: OTI1F3]
Reactive Species	Dog, Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-OXSR1 mouse monoclonal antibody, clone OTI1F3 (formerly 1F3). Catalog# M05141. Tested in IHC, WB. This antibody reacts with Human, Monkey, Mouse, Rat, Dog.
Application	IHC, WB
Clonality	Monoclonal OTI1F3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O95747


### Technical Details

Immunogen	Full-length protein expressed in 293T cell transfected with human OXSR1 expression vector
Isotype	IgG1
Concentration	1.28 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:1000 IHC 1:150


## Anti-OXSR1 Mouse Monoclonal Antibody [Clone ID: OTI1F3] (M05141) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY OXSR1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-OXSR1.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-OXSR1 monoclonal antibody.



Immunohistochemical staining of paraffin-embedded Human liver tissue within the normal limits using anti-OXSR1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-OXSR1 Mouse Monoclonal Antibody [Clone ID: OTI1F3]

For Research Use Only. Not for use in diagnostic procedures.