

## Anti-ULBP3 Antibody (Monoclonal, 33U90)

Catalog Number: M05191

### About ULBP3

ULBP3 (UL16 binding protein-3), a class 1 major histocompatibility complex-like (MHC-like) molecule and member of the family of ligands of the killer cell lectin-like receptor subfamily K member 1 (NKG2D) (PMID:31292299). ULBP3 are expressed at the cell surface of multiple types of cells and tissues as a glycosylphosphatidylinositol (GPI)-linked proteins (PMID:22741585).

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-ULBP3 Antibody (Monoclonal, 33U90)   |
| Reactive Species     | Human   |
| Description          | Boster Bio Anti-ULBP3 Antibody (Monoclonal, 33U90) catalog # M05191. Tested in WB, IHC, IP applications. This antibody reacts with Human.   |
| Application          | IP, IHC, WB   |
| Clonality            | Monoclonal 33U90  |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.  |
| Host                 | Rabbit  |
| Uniprot ID           | Q9BZM4  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Recombinant protein within human ULBP3 aa 23-218.                                       |
| Form                | Liquid  |
| Concentration       | 500 ug/ml   |
| Purification        | Protein A affinity purified.  |
| Suggested Dilutions | Western blot, 1:500-2000<br>Immunohistochemistry, 1:50-200<br>ImmunoPrecipitation, 1:50 |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ULBP3 Antibody (Monoclonal, 33U90)

For Research Use Only. Not for use in diagnostic procedures.