

Anti-HDAC10 Rabbit Monoclonal Antibody

Catalog Number: M05215


Overview

Product Name	Anti-HDAC10 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-HDAC10 Rabbit Monoclonal Antibody catalog # M05215. Tested in WB, ICC/IF, IP applications. This antibody reacts with Human, Mouse.
Application	IP, IF, ICC, WB
Clonality	Monoclonal HEO-8
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q969S8


Technical Details

Immunogen	A synthesized peptide derived from human HDAC10
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:1000-1:2000 ICC/IF 1:50-1:200 IP 1:50


Anti-HDAC10 Rabbit Monoclonal Antibody (M05215) Images



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of HDAC10 expression in (1) HeLa cell lysate; (2) 3T3 cell lysate.



Immunofluorescent analysis of HeLa cells, using HDAC10 Antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HDAC10 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.