

Anti-HDAC10 Mouse Monoclonal Antibody [Clone ID: OTI8G5]

Catalog Number: M05215-1


Overview

Product Name	Anti-HDAC10 Mouse Monoclonal Antibody [Clone ID: OTI8G5]
Reactive Species	Human
Description	Boster Bio Anti-HDAC10 mouse monoclonal antibody, clone OTI8G5 (formerly 8G5). Catalog# M05215-1. Tested in IF, IHC, WB. This antibody reacts with Human.
Application	IF, IHC, WB
Clonality	Monoclonal OTI8G5
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q969S8


Technical Details

Immunogen	Full length human recombinant protein of human HDAC10 (NP_114408) produced in HEK293T cell.
Isotype	IgG1
Concentration	0.76 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:1000 IHC 1:50 IF 1:100


Anti-HDAC10 Mouse Monoclonal Antibody [Clone ID: OTI8G5] (M05215-1) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HDAC10 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HDAC10.



Immunohistochemical staining of paraffin-embedded Carcinoma of prostate tissue using anti-HDAC10 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)





Immunohistochemical staining of paraffin-embedded Adenocarcinoma of colon tissue using anti-HDAC10 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)




Immunohistochemical staining of paraffin-embedded Kidney tissue within the normal limits using anti-HDAC10 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)

Immunohistochemical staining of paraffin-embedded liver tissue within the normal limits using anti-HDAC10 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Carcinoma of liver tissue using anti-HDAC10 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Anti-HDAC10 mouse monoclonal antibody (M05215-1) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HDAC10.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HDAC10 Mouse Monoclonal Antibody [Clone ID: OTI8G5]

For Research Use Only. Not for use in diagnostic procedures.