

## Anti-Cbx8 Antibody Picoband® (monoclonal, 8G7) Cy3 Conjugated

Catalog Number: M05234-Cy3

### About CBX8

CBX8 functions as a transcriptional repressor and has a role in DNA damage response. This gene is mapped to chromosome 17q25.3 based on an alignment of the CBX8 sequence with the genomic sequence (GRCh38).

### Overview

Product Name	Anti-Cbx8 Antibody Picoband® (monoclonal, 8G7) Cy3 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Monoclonal 8G7
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	Q9HC52

### Technical Details

Immunogen	E.coli-derived human Cbx8 recombinant protein (Position: M1-R389).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2b
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cbx8 Antibody (monoclonal, 8G7) - Cy3

For Research Use Only. Not for use in diagnostic procedures.