

Anti-TSPAN26 / CD37 Reference Antibody (AGS67E)

Catalog Number: M05276-1


Overview

| | |
|----------------------|--|
| Product Name | Anti-TSPAN26 / CD37 Reference Antibody (AGS67E) |
| Reactive Species | Cynomolgus monkey, Human |
| Description | Boster Bio Anti-TSPAN26 / CD37 Reference Antibody (AGS67E) (Catalog # M05276-1). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.880EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application | ELISA, Flow Cytometry, Functional Assay, Kinetics |
| Clonality | Monoclonal |
| Formulation | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles. |
| Uniprot ID | P11049 |


Technical Details

| | |
|---------|--------|
| Isotype | IgG2SA |
| Form | Liquid |

Anti-TSPAN26 / CD37 Reference Antibody (AGS67E) (M05276-1) Images



Anti-TSPAN26/CD37 Reference Antibody (AGS67E) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TSPAN26/CD37 Reference Antibody (AGS67E) is more than 100%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TSPAN26 / CD37 Reference Antibody (AGS67E)

For Research Use Only. Not for use in diagnostic procedures.