

## Anti-TSPAN26 / CD37 Reference Antibody (lilotomab)

Catalog Number: M05276-2


### Overview

Product Name	Anti-TSPAN26 / CD37 Reference Antibody (lilotomab)
Reactive Species	Human
Description	Boster Bio Anti-TSPAN26 / CD37 Reference Antibody (lilotomab) (Catalog # M05276-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.881EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P11049


### Technical Details

Isotype	IgG1
Form	Liquid

## Anti-TSPAN26 / CD37 Reference Antibody (lilotomab) (M05276-2) Images



Anti-TSPAN26/CD37 Reference Antibody (lilotomab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TSPAN26/CD37 Reference Antibody (lilotomab) is more than 100%

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TSPAN26 / CD37 Reference Antibody (lilotomab)

For Research Use Only. Not for use in diagnostic procedures.