

Anti-(DANRE) srrt Antibody (C-Term)

Catalog Number: M05444

About srrt

Acts as a mediator between the cap-binding complex (CBC) and the primary microRNAs (miRNAs) processing machinery during cell proliferation. Contributes to the stability and delivery of capped primary miRNA transcripts to the primary miRNA processing complex, thereby playing a role in RNA-mediated gene silencing (RNAi) by miRNAs (By similarity). Independently of its activity on miRNAs, necessary and sufficient to promote neural stem cell self-renewal. Does so by directly binding SOX2 promoter and positively regulating its transcription (By similarity).


Overview

Product Name	Anti-(DANRE) srrt Antibody (C-Term)
Reactive Species	Human, Zebrafish
Description	Boster Bio Anti-(DANRE) srrt Antibody (C-Term) (Catalog # M05444). Tested in WB application(s). This antibody reacts with Human, Zebrafish.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q66I22

Technical Details

Immunogen	This DANRE srrt antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 733-766 amino acids from DANRE srrt.
Predicted Reactive Species	Rat
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:2000

Anti-(DANRE) srvt Antibody (C-Term) (M05444) Images



All lanes : Anti-srvt Antibody (C-Term) at 1:2000 dilution
Lane 1: Zebrafish muscle lysates
Lane 2: Zebrafish lysates
Lane 3: Jurkat whole cell lysates
Lysates/proteins at 20 µg per lane.
Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 103 kDa
Blocking/Dilution buffer: 5% NFDN/TBST.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-(DANRE) srvt Antibody (C-Term)

For Research Use Only. Not for use in diagnostic procedures.