

Anti-DBN1 Antibody (Center)

Catalog Number: M05530-2

About DBN1

Drebrins might play some role in cell migration, extension of neuronal processes and plasticity of dendrites. Required for actin polymerization at immunological synapses (IS) and for CXCR4 recruitment to IS.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-DBN1 Antibody (Center)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-DBN1 Antibody (Center) (Catalog # M05530-2). Tested in WB application(s). This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16643


Technical Details

Immunogen	This DBN1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 194-228 amino acids from the Central region of human DBN1.
Predicted Reactive Species	Rat
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.

Suggested Dilutions

WB: 1:2000

Anti-DBN1 Antibody (Center) (M05530-2) Images



All lanes : Anti-DBN1 Antibody (Center) at 1:2000 dilution
Lane 1: HeLa whole cell lysate
Lane 2: Jurkat whole cell lysate
Lane 3: SH-SY5Y whole cell lysate
Lane 4: mouse brain lysate
Lane 5: Neuro-2a whole cell lysate
Lysates/proteins at 20 µg per lane. Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 71 kDa
Blocking/Dilution buffer: 5% NFDM/TBST.

Anti-DBN1 Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.