

## Anti-Drebrin/DBN1 Antibody Picoband® (monoclonal, 4F6E7) APC Conjugated

Catalog Number: M05530-4-APC

### About DBN1

Drebrin is a protein that in humans is encoded by the DBN1 gene. The protein encoded by this gene is a cytoplasmic actin-binding protein thought to play a role in the process of neuronal growth. It is a member of the drebrin family of proteins that are developmentally regulated in the brain. A decrease in the amount of this protein in the brain has been implicated as a possible contributing factor in the pathogenesis of memory disturbance in Alzheimer's disease. At least two alternative splice variants encoding different protein isoforms have been described for this gene.

### Overview

Product Name	Anti-Drebrin/DBN1 Antibody Picoband® (monoclonal, 4F6E7) APC Conjugated
Reactive Species	Human, Monkey, Mouse, Rat
Application	Flow Cytometry
Clonality	Monoclonal 4F6E7
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	Q16643

### Technical Details

Immunogen	E.coli-derived human Drebrin/DBN1 recombinant protein (Position: H9-D649).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2a
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Drebrin/DBN1 Antibody (monoclonal, 4F6E7) - APC

For Research Use Only. Not for use in diagnostic procedures.