

## Anti-PPCS Antibody Picoband® (monoclonal, 7G13) Fluoro594 Conjugated

Catalog Number: M05567-Fluoro594

### About PPCS

Phosphopantothenoylcysteine synthetase in humans is encoded by the PPCS gene. Biosynthesis of coenzyme A (CoA) from pantothenic acid (vitamin B5) is an essential universal pathway in prokaryotes and eukaryotes. PPCS, one of the last enzymes in this pathway, converts phosphopantothenate to phosphopantothenoylcysteine. By genomic sequence analysis, this PPCS gene is mapped to chromosome 1.

### Overview

Product Name	Anti-PPCS Antibody Picoband® (monoclonal, 7G13) Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Monoclonal 7G13
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	Q9HAB8

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human PPCS, which shares 73.9% and 78.2% amino acid (aa) sequence identity with mouse and rat PPCS, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2a
Form	Liquid
Concentration	0.5 mg/mL
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PPCS Antibody (monoclonal, 7G13) - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.