

Anti-Nidogen 1 Rabbit Monoclonal Antibody

Catalog Number: M05621


Overview

Product Name	Anti-Nidogen 1 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Nidogen 1 Rabbit Monoclonal Antibody catalog # M05621. Tested in WB application. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 25N81
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P14543

Technical Details

Immunogen	A synthesized peptide derived from human Nidogen 1
Isotype	IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000

Anti-Nidogen 1 Rabbit Monoclonal Antibody (M05621) Images



Western blot analysis of Nidogen 1 expression in Human fetal heart lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nidogen 1 Rabbit Monoclonal Antibody