

Anti-Napsin-A NAPSA Rabbit Monoclonal Antibody, Clone#RM366

Catalog Number: M05685-3


Overview

Product Name	Anti-Napsin-A NAPSA Rabbit Monoclonal Antibody, Clone#RM366
Reactive Species	Human
Description	Boster Bio Anti-Napsin-A NAPSA Rabbit Monoclonal Antibody, Clone#RM366 (Catalog # M05685-3). Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal RM366
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O96009


Technical Details

Immunogen	A peptide corresponding to the N-terminus of human Napsin-A
Cross Reactivity	This antibody reacts human Napsin-A
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB: 1:1000-1:2000 dilution

Anti-Napsin-A NAPSA Rabbit Monoclonal Antibody, Clone#RM366 (M05685-3) Images



Western Blotting resultWestern Blot of human lung tissue lysate using anti-Napsin-A rabbit monoclonal antibody (Clone RM366) at a 1:2000 dilution.



Western Blotting resultWestern Blot of human lung tissue lysate using anti-Napsin-A rabbit monoclonal antibody (Clone RM366) at a 1:2000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Napsin-A NAPSA Rabbit Monoclonal Antibody, Clone#RM366

For Research Use Only. Not for use in diagnostic procedures.