

## Anti-SCP-3 SYCP3 Monoclonal Antibody

Catalog Number: M05718


### Overview

Product Name	Anti-SCP-3 SYCP3 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-SCP-3 SYCP3 Monoclonal Antibody catalog # M05718. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 3B10
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q8IZU3

### Technical Details

Immunogen	Purified recombinant human SCP-3 protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000

## Anti-SCP-3 SYCP3 Monoclonal Antibody (M05718) Images



Western Blot analysis using SCP-3 Monoclonal Antibody against mouse pancreas lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SCP-3 SYCP3 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.