

## Anti-DUSP5 Mouse Monoclonal Antibody [Clone ID: OTI5H9]

Catalog Number: M05809


### Overview

Product Name	Anti-DUSP5 Mouse Monoclonal Antibody [Clone ID: OTI5H9]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio DUSP5 mouse monoclonal antibody, clone OTI5H9. Catalog# M05809. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI5H9
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q16690


### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 125-384 of human DUSP5 (NP_004410) produced in E.coli.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500~2000

## Anti-DUSP5 Mouse Monoclonal Antibody [Clone ID: OTI5H9] (M05809) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DUSP5 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DUSP5 (1:2000).



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-DUSP5 monoclonal antibody (1:500).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DUSP5 Mouse Monoclonal Antibody [Clone ID: OTI5H9]

For Research Use Only. Not for use in diagnostic procedures.