

Anti-UBR2 Picoband® Antibody (monoclonal, 2G10) Biotin Conjugated

Catalog Number: M05812-Biotin

About UBR2

E3 ubiquitin-protein ligase UBR2 is an enzyme that in humans is encoded by the UBR2 gene. It is mapped to 6p21.1. This gene encodes an E3 ubiquitin ligase of the N-end rule proteolytic pathway that targets proteins with destabilizing N-terminal residues for polyubiquitylation and proteasome-mediated degradation. Alternative splicing results in multiple transcript variants.

Overview

Product Name	Anti-UBR2 Picoband® Antibody (monoclonal, 2G10) Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Monoclonal 2G10
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Mouse
Uniprot ID	Q8IWW8

Technical Details

Immunogen	E.coli-derived human UBR2 recombinant protein (Position: Q10-A664).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2a
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBR2 Antibody (monoclonal, 2G10) - Biotin

For Research Use Only. Not for use in diagnostic procedures.