

## Anti-UBR2 Picoband® Antibody (monoclonal, 2G10) Fluoro488 Conjugated

Catalog Number: M05812-Fluoro488

### About UBR2

E3 ubiquitin-protein ligase UBR2 is an enzyme that in humans is encoded by the UBR2 gene. It is mapped to 6p21.1. This gene encodes an E3 ubiquitin ligase of the N-end rule proteolytic pathway that targets proteins with destabilizing N-terminal residues for polyubiquitylation and proteasome-mediated degradation. Alternative splicing results in multiple transcript variants.

### Overview

Product Name	Anti-UBR2 Picoband® Antibody (monoclonal, 2G10) Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Monoclonal 2G10
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	Q8IWW8

### Technical Details

Immunogen	E.coli-derived human UBR2 recombinant protein (Position: Q10-A664).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2a
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBR2 Antibody (monoclonal, 2G10) - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.