

## Anti-PGLYRP1 / PGRP-S Reference Antibody (Novo Nordisk patent anti-PGLYRP1)

Catalog Number: M05914


### Overview

Product Name	Anti-PGLYRP1 / PGRP-S Reference Antibody (Novo Nordisk patent anti-PGLYRP1)
Reactive Species	Human
Description	Boster Bio Anti-PGLYRP1 / PGRP-S Reference Antibody (Novo Nordisk patent anti-PGLYRP1) (Catalog # M05914). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.898EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	O75594


### Technical Details

Isotype	IgG4SP
Form	Liquid

## Anti-PGLYRP1 / PGRP-S Reference Antibody (Novo Nordisk patent anti-PGLYRP1) (M05914) Images



Anti-PGLYRP1/PGRP-S Reference Antibody (Novo Nordisk patent anti-PGLYRP1) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-PGLYRP1/PGRP-S Reference Antibody (Novo Nordisk patent anti-PGLYRP1) is more than 95%

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PGLYRP1 / PGRP-S Reference Antibody (Novo Nordisk patent anti-PGLYRP1)

For Research Use Only. Not for use in diagnostic procedures.