

## Anti-Carboxypeptidase A1+A2+B Rabbit Monoclonal Antibody

Catalog Number: M05985-2


### Overview

Product Name	Anti-Carboxypeptidase A1+A2+B Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Carboxypeptidase A1+A2+B Rabbit Monoclonal Antibody catalog # M05985-2. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Monoclonal 23C13
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P15085/P15086/P48052

### Technical Details

Immunogen	A synthesized peptide derived from human Carboxypeptidase A1+A2+B
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200

## Anti-Carboxypeptidase A1+A2+B Rabbit Monoclonal Antibody (M05985-2) Images



IHC analysis of Carboxy peptidase A1/A2/B1 using anti-Carboxy peptidase A1/A2/B1 antibody (M05985-2). Carboxy peptidase A1/A2/B1 was detected in a paraffin-embedded section of human pancreas tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 1:50 rabbit anti-Carboxy peptidase A1/A2/B1 Antibody (M05985-2) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Carboxypeptidase A1+A2+B Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.