

Anti-Rab2A (C-Terminus) monoclonal Antibody

Catalog Number: M06007


Overview

Product Name	Anti-Rab2A (C-Terminus) monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Rab2A (C-Terminus) monoclonal Antibody catalog # M06007. Tested in WB,ICC applications. This antibody reacts with Human.
Application	ICC, WB
Clonality	Monoclonal 3E9
Formulation	Mouse IgG1, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P61019


Technical Details

Immunogen	Synthetic peptide, corresponding to the C-Terminus of Human Rab2A.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:800 ICC: 1:100-200

Anti-Rab2A (C-Terminus) monoclonal Antibody (M06007) Images



Western blot detection of RAB2A in HeLa, Ramos&K562 cell lysates and using RAB2A antibody (1:800 diluted) .



Immunocytochemistry of HeLa cells using anti- RAB2A(C-Terminus) antibody diluted 1:200

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rab2A (C-Terminus) monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.