

Anti-Rab 2A Monoclonal Antibody

Catalog Number: M06007-1


Overview

Product Name	Anti-Rab 2A Monoclonal Antibody
Reactive Species	Bovine, Chicken, Dog, Human, Mouse, Pig, Rat, Zebrafish
Description	Boster Bio Anti-Rab 2A Monoclonal Antibody catalog # M06007-1. Tested in IF, WB applications. This antibody reacts with Bovine, Canine, Chicken, Human, Mouse, Rat, Swine, Zebrafish.
Application	IF, WB
Clonality	Monoclonal 3E9
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P61019


Technical Details

Immunogen	Purified recombinant human Rab 2A (C-terminus) protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000 IF 1:100-1:500

Anti-Rab 2A Monoclonal Antibody (M06007-1) Images



Western Blot analysis using Rab 2A Monoclonal Antibody against HeLa cell lysate.



Immunofluorescence analysis of HeLa cells using Rab 2A Monoclonal Antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rab 2A Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.