

## Anti-TNFSF9 / 4-1BBL Reference Antibody (Abbvie patent anti-TNFSF9)

Catalog Number: M06032-2

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




### Overview

|                      |                                                                                                                                                                                                                                                           |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-TNFSF9 / 4-1BBL Reference Antibody (Abbvie patent anti-TNFSF9)                                                                                                                                                                                       |
| Reactive Species     | Human                                                                                                                                                                                                                                                     |
| Description          | Boster Bio Anti-TNFSF9 / 4-1BBL Reference Antibody (Abbvie patent anti-TNFSF9) (Catalog # M06032-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.900EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application          | ELISA, Flow Cytometry, Functional Assay, Kinetics                                                                                                                                                                                                         |
| Clonality            | Monoclonal                                                                                                                                                                                                                                                |
| Formulation          | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.                                                                                                                           |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.                                                                                                                 |
| Uniprot ID           | P41273                                                                                                                                                                                                                                                    |


### Technical Details

|         |        |
|---------|--------|
| Isotype | IgG1   |
| Form    | Liquid |

## Anti-TNFSF9 / 4-1BBL Reference Antibody (Abbvie patent anti-TNFSF9) (M06032-2) Images



Anti-TNFSF9/4-1BBL Reference Antibody (Abbvie patent anti-TNFSF9) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-TNFSF9/4-1BBL Reference Antibody (Abbvie patent anti-TNFSF9) is more than 95%

### Anti-TNFSF9 / 4-1BBL Reference Antibody (Abbvie patent anti-TNFSF9)

For Research Use Only. Not for use in diagnostic procedures.