

Anti-EIF3J Antibody (Monoclonal, 32E12)

Catalog Number: M06040

About EIF3J

EIF3 has a key role in binding of initiator methionyl-tRNA and mRNA to the 40S ribosomal subunit to form the 40S initiation complex (PMID:17588516). The eIF3 complex stimulates several steps in the translation initiation pathway, including dissociation of 80S ribosomes into 40S and 60S subunits, binding of a ternary complex (TC) consisting of Met-tRNA, eIF2, and GTP to the small subunit (forming the 43S preinitiation complex) and recruitment of mRNA to the 43S complex to produce the 48S complex (PMID:11560931). EIF3J was identified as a high copy suppressor of the temperature-sensitive (Ts-) phenotype of the rpg1-1 allele of TIF32/RPG1, encoding the largest subunit of yeast eIF3 (eIF3a) (PMID:10488093).

Overview

Product Name	Anti-EIF3J Antibody (Monoclonal, 32E12)
Reactive Species	Human, Monkey, Rat
Description	Boster Bio Anti-EIF3J Antibody (Monoclonal, 32E12) catalog # M06040. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human, Rat, Monkey.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal 32E12
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75822

Technical Details

Immunogen	Recombinant protein within human EIF3J aa 3-256.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.

Suggested Dilutions

Western blot, 1:500-2000
Immunohistochemistry, 1:50-200
Immunocytochemistry/Immunofluorescence, 1:50-200
Flow Cytometry (Fixed), 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EIF3J Antibody (Monoclonal, 32E12)

For Research Use Only. Not for use in diagnostic procedures.