

Anti-NIRF (UHRF2) Mouse Monoclonal Antibody [Clone ID: OTI4E10]

Catalog Number: M06294


Overview

Product Name	Anti-NIRF (UHRF2) Mouse Monoclonal Antibody [Clone ID: OTI4E10]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio UHRF2 mouse monoclonal antibody, clone OTI4E10. Catalog# M06294. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI4E10
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q96PU4


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 1-288 of human UHRF2 (NP_690856) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB: 1:500~2000

Anti-NIRF (UHRF2) Mouse Monoclonal Antibody [Clone ID: OTI4E10] (M06294) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UHRF2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UHRF2 (1:2000).



Western blot analysis of extracts (35ug) from 7 different cell lines by using anti-UHRF2 monoclonal antibody (1:500).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NIRF (UHRF2) Mouse Monoclonal Antibody [Clone ID: OTI4E10]

For Research Use Only. Not for use in diagnostic procedures.