

Anti-LZTS1 Mouse Monoclonal Antibody [Clone ID: OTI3D12]

Catalog Number: M06359


Overview

Product Name	Anti-LZTS1 Mouse Monoclonal Antibody [Clone ID: OTI3D12]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio LZTS1 mouse monoclonal antibody, clone OTI3D12. Catalog# M06359. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI3D12
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9Y250


Technical Details

Immunogen	Full length human recombinant protein of human LZTS1 (NP_066300) produced in HEK293T Cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500~2000

Anti-LZTS1 Mouse Monoclonal Antibody [Clone ID: OTI3D12] (M06359) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY LZTS1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-LZTS1 (1:2000).



Western blot analysis of extracts (35ug) from Jurkat cell lines by using anti-LZTS1 monoclonal antibody (1:500).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LZTS1 Mouse Monoclonal Antibody [Clone ID: OTI3D12]

For Research Use Only. Not for use in diagnostic procedures.