

Anti-LI-Cadherin-17 CDH17-Rabbit Monoclonal Antibody

Catalog Number: M06389


Overview

Product Name	Anti-LI-Cadherin-17 CDH17-Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-LI-Cadherin-17 CDH17-Rabbit Monoclonal Antibody catalog # M06389. Tested in WB, IHC, IP applications. This antibody reacts with Human.
Application	IP, IHC, WB
Clonality	Monoclonal AFFC-3
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q12864


Technical Details

Immunogen	A synthesized peptide derived from human LI-Cadherin
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 IP 1:50

Anti-LI-Cadherin-17 CDH17-Rabbit Monoclonal Antibody (M06389) Images



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of LI-Cadherin expression in HT-29 cell lysate.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LI-Cadherin-17 CDH17-Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.