

## Anti-VPS28 Mouse Monoclonal Antibody [Clone ID: OTI1A8]

Catalog Number: M06407


### Overview

Product Name	Anti-VPS28 Mouse Monoclonal Antibody [Clone ID: OTI1A8]
Reactive Species	Dog, Human, Monkey, Mouse, Rat
Description	Boster Bio VPS28 mouse monoclonal antibody, clone OTI1A8 (formerly 1A8). Catalog# M06407. Tested in IF, IHC, WB. This antibody reacts with Human, Monkey, Mouse, Rat, Dog.
Application	IF, IHC, WB
Clonality	Monoclonal OTI1A8
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9UK41


### Technical Details

Immunogen	Full length human recombinant protein of human VPS28 (NP_057292) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB: 1:200~2000 IHC: 1:150 IF: 1:100


## Anti-VPS28 Mouse Monoclonal Antibody [Clone ID: OTI1A8] (M06407) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY VPS28 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-VPS28.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-VPS28 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



Western blot analysis of extracts (10ug) from a mouse cell line and 3 different mouse tissues by using anti-VPS28 monoclonal antibody (1:200).



Immunohistochemical staining of paraffin-embedded Human pancreas tissue within the normal limits using anti-VPS28 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Anti-VPS28 mouse monoclonal antibody (M06407) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY VPS28.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-VPS28 Mouse Monoclonal Antibody [Clone ID: OT11A8]

For Research Use Only. Not for use in diagnostic procedures.