

Anti-FGFBP1 Mouse Monoclonal Antibody [Clone ID: OTI2B6]

Catalog Number: M06468


Overview

Product Name	Anti-FGFBP1 Mouse Monoclonal Antibody [Clone ID: OTI2B6]
Reactive Species	Human
Description	Boster Bio FGFBP1 mouse monoclonal antibody, clone OTI2B6 (formerly 2B6). Catalog# M06468. Tested in IHC, WB. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal OTI2B6
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q14512


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 24-234 of human FGFBP1 (NP_005121) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IHC 1:150


Anti-FGFBP1 Mouse Monoclonal Antibody [Clone ID: OTI2B6] (M06468) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FGFBP1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FGFBP1 (1:2000).



Immunohistochemical staining of paraffin-embedded Human Ovary tissue within the normal limits using anti-FGFBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by Tris-EDTA)



Immunohistochemical staining of paraffin-embedded Human prostate tissue within the normal limits using anti-FGFBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by Tris-EDTA)

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGFBP1 Mouse Monoclonal Antibody [Clone ID: OTI2B6]

For Research Use Only. Not for use in diagnostic procedures.