

Anti-PDLIM2 Mouse Monoclonal Antibody [Clone ID: OTI11A12]

Catalog Number: M06506


Overview

Product Name	Anti-PDLIM2 Mouse Monoclonal Antibody [Clone ID: OTI11A12]
Reactive Species	Human
Description	Boster Bio PDLIM2 mouse monoclonal antibody, clone OTI11A12 (formerly 11A12). Catalog# M06506. Tested in FC, IF, IHC, WB. This antibody reacts with Human.
Application	Flow Cytometry, IF, IHC, WB
Clonality	Monoclonal OTI11A12
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q96JY6


Technical Details

Immunogen	Full length human recombinant protein of human PDLIM2 (NP_067643) produced in HEK293T cell.
Isotype	IgG1
Concentration	0.62 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IHC 1:100 IF 1:100 Flow Cytometry 1:100


Anti-PDLIM2 Mouse Monoclonal Antibody [Clone ID: OTI11A12] (M06506) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PDLIM2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PDLIM2.




Immunohistochemical staining of paraffin-embedded Human lymph node tissue within the normal limits using anti-PDLIM2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)




Immunohistochemical staining of paraffin-embedded Human lymphoma tissue using anti-PDLIM2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)




Immunohistochemical staining of paraffin-embedded Human Kidney tissue within the normal limits using anti-PDLIM2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)




Immunohistochemical staining of paraffin-embedded Human liver tissue within the normal limits using anti-PDLIM2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)




HEK293T cells transfected with either PDLIM2 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-PDLIM2 antibody (M06506)



Flow cytometric Analysis of HeLa cells



Flow cytometric Analysis of Jurkat cells



Anti-PDLIM2 mouse monoclonal antibody (M06506) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PDLIM2.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PDLIM2 Mouse Monoclonal Antibody [Clone ID: OT111A12]

For Research Use Only. Not for use in diagnostic procedures.