

Anti-Uroplakin3a UPK3A Monoclonal Antibody

Catalog Number: M06534-1


Overview

Product Name	Anti-Uroplakin3a UPK3A Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Uroplakin3a UPK3A Monoclonal Antibody catalog # M06534-1. Tested in IHC, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IHC
Clonality	Monoclonal AFEH-21
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75631


Technical Details

Immunogen	A synthesized peptide derived from human Uroplakin III Component of the asymmetric unit membrane (AUM) ; a highly specialized biomembrane elaborated by terminally differentiated urothelial cells. May play an important role in AUM-cytoskeleton interaction in terminally differentiated urothelial cells.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	IHC 1:1000-1:500 FC 1:100


Anti-Uroplakin3a UPK3A Monoclonal Antibody (M06534-1) Images




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Immunohistochemical analysis of paraffin-embedded human bladder, using Uroplakin III Antibody.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Uroplakin3a UPK3A Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.