

Anti-EIF4A2 Antibody (Monoclonal, 33E77)

Catalog Number: M06549-1

About EIF4A2

Eukaryotic initiation factor 4A-II is a protein that in humans is encoded by the EIF4A2 gene. It is mapped to 18p11.2. Eukaryotic initiation factor 4A plays an important role in the binding of mRNA to the 43S preinitiation complex when protein synthesis begins. Two highly homologous forms of functional EIF4A genes, Eif4a1 and Eif4a2, have been isolated in mice; yeast cells also possess 2 EIF4A genes, TIF1 and TIF2. The murine Eif4a and yeast TIF genes appear to belong to a DEAD-box gene family, whose members exhibit extensive amino acid similarity and contain the asp-glu-ala-asp (DEAD) sequence. DEAD-box genes have been identified in species ranging from E-coli to humans. Their function appears to be related to transcriptional/translational regulation.

Overview

Product Name	Anti-EIF4A2 Antibody (Monoclonal, 33E77)
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-EIF4A2 Antibody (Monoclonal, 33E77) catalog # M06549-1. Tested in WB, IP applications. This antibody reacts with Human, Mouse, Rat, Monkey.
Application	IP, WB
Clonality	Monoclonal 33E77
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14240

Technical Details

Immunogen	Recombinant protein within Human eIF4A2 aa 3-405.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.

Suggested Dilutions

Western blot, 1:500-2000
ImmunoPrecipitation, 1:50

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EIF4A2 Antibody (Monoclonal, 33E77)

For Research Use Only. Not for use in diagnostic procedures.