

## Anti-Semaphorin 3B/SEMA3B Antibody Picoband® (monoclonal, 9C4F7) Fluoro594 Conjugated

Catalog Number: M06559-Fluoro594

### About SEMA3B

Semaphorin-3B is a protein that in humans is encoded by the SEMA3B gene. It is mapped to 3p21.31. The semaphorin/collapsin family of molecules plays a critical role in the guidance of growth cones during neuronal development. The secreted protein encoded by this gene family member is important in axonal guidance and has been shown to act as a tumor suppressor by inducing apoptosis.

### Overview

Product Name	Anti-Semaphorin 3B/SEMA3B Antibody Picoband® (monoclonal, 9C4F7) Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Monoclonal 9C4F7
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	Q13214

### Technical Details

Immunogen	E.coli-derived human Semaphorin 3B/SEMA3B recombinant protein (Position: Y55-R624).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2b
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Semaphorin 3B/SEMA3B Antibody (monoclonal, 9C4F7) - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.