

## Anti-RAB5B Antibody

Catalog Number: M06608

### About RAB5B

Protein transport. Probably involved in vesicular traffic (By similarity).


### Overview

Product Name	Anti-RAB5B Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-RAB5B Antibody (Catalog # M06608). Tested in WB application(s). This antibody reacts with Human, Mouse.
Application	WB
Clonality	Monoclonal 1615CT175.80.21
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	P61020


### Technical Details

Immunogen	This RAB5B antibody is generated from a mouse immunized with a recombinant protein of human RAB5B.
Predicted Reactive Species	Mouse
Isotype	IgG1,k
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	WB: 1:2000

## Anti-RAB5B Antibody (M06608) Images



Anti-RAB5B Antibody at 1:2000 dilution + A431 whole cell lysateLysates/proteins at 20 µg per lane. SecondaryGoat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 24 kDaBlocking/Dilution buffer: 5% NFDM/TBST.



All lanes : Anti-RAB5B Antibody at 1:2000 dilutionLane 1: HeLa whole cell lysateLane 2: HepG2 whole cell lysateLane 3: Jurkat whole cell lysateLane 4: NIH-3T3 whole cell lysateLysates/proteins at 20 µg per lane. SecondaryGoat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 24 kDaBlocking/Dilution buffer: 5% NFDM/TBST.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAB5B Antibody

For Research Use Only. Not for use in diagnostic procedures.