

Anti-PCSK7 Mouse Monoclonal Antibody [Clone ID: OTI1B8]

Catalog Number: M06686


Overview

Product Name	Anti-PCSK7 Mouse Monoclonal Antibody [Clone ID: OTI1B8]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio PCSK7 mouse monoclonal antibody, clone OTI1B8 (formerly 1B8). Catalog# M06686. Tested in FC, IF. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF
Clonality	Monoclonal OTI1B8
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q16549


Technical Details

Immunogen	Full length human recombinant protein of human PCSK7 (NP_004707) produced in HEK293T cell.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	IF 1:100 Flow Cytometry 1:100

Anti-PCSK7 Mouse Monoclonal Antibody [Clone ID: OTI1B8] (M06686) Images



HEK293T cells transfected with either PCSK7 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-PCSK7 antibody (M06686)



Anti-PCSK7 mouse monoclonal antibody (M06686) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PCSK7.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PCSK7 Mouse Monoclonal Antibody [Clone ID: OTI1B8]

For Research Use Only. Not for use in diagnostic procedures.