

Anti-TRMT1 (TRMU) Mouse Monoclonal Antibody [Clone ID: OTI2A4]

Catalog Number: M06698


Overview

Product Name	Anti-TRMT1 (TRMU) Mouse Monoclonal Antibody [Clone ID: OTI2A4]
Reactive Species	Dog, Human, Monkey, Mouse, Rat
Description	Boster Bio TRMU mouse monoclonal antibody, clone OTI2A4 (formerly 2A4). Catalog# M06698. Tested in IF, WB. This antibody reacts with Human, Monkey, Mouse, Rat, Dog.
Application	IF, WB
Clonality	Monoclonal OTI2A4
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O75648


Technical Details

Immunogen	Full length human recombinant protein of human TRMU (NP_060476) produced in HEK293T cell.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/


Anti-TRMT1 (TRMU) Mouse Monoclonal Antibody [Clone ID: OTI2A4] (M06698) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TRMU (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRMU.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-TRMU monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



Anti-TRMU mouse monoclonal antibody (M06698) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TRMU.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRMT1 (TRMU) Mouse Monoclonal Antibody [Clone ID: OT12A4]